



Tough Enough.

We passed all the tests. No one else has.

Engineered to perform in adverse and demanding conditions.
The iCRco 3600 (ClearVet CR) is the US Military's system of choice.

Independently Tested and Certified for Environmental Durability and Ruggedness

Some of the tests the iCRco 3600 (ClearVet CR) passed and others failed:

- VIBRATION TESTING—10 minutes per axis. PASSED
- SHOCK TESTING—30G Shock—18 Repetitions. PASSED
- TRANSIT DROP TEST—10 drops from two feet. PASSED
- BLOWING SAND TEST—Blowing Sand at 45 MPH for hours. PASSED
- SETTLING DUST—6 gm/m²/day for three days. PASSED
- HUMIDITY TEST—5 day extreme humidity cycle. PASSED
- SALT FOG TEST—5% Sodium Chloride Fog for 48 hours. PASSED
- FUNGUS TEST—28 day exposure to Fungus Spores. PASSED
- TEMPERATURE TESTING—Temperatures from of 30 degrees F to 105 degrees F. PASSED

The system remained operational through ALL tests

**Durability makes the iCRco 3600 (ClearVet CR) the superior system for your
Mobile Veterinary application**

