

# apexx AE-1020HF High Frequency Portable X-Ray

The apexx AE-1020HF High Frequency Portable X-Ray is an efficient high frequency portable x-ray that is light weight, durable and powerful.

The small x-ray tube focal spot produces an excellent image with fine detail and superb resolution.



AC Input .....	VAC - 110VAC, 50/60 Hz
X-ray Tube .....	Toshiba
Focal Spot Size .....	1.2 mm
Anode Heat Storage .....	20 kHU
Total Filtration .....	2.6 mm Al
Collimator .....	Dual knob, fully adjustable
Exposure Switch.....	Two stage
Exposure Cord .....	9 feet
Power Cord .....	15 feet
Size.....	9"W x 7"H x 12"L
Weight.....	24.2 lb (11 kg)

Output:

kVp Range (2kV steps)	mAs Range
40 – 48 kVp	0.4 – 50 mAs
50 – 58 kVp	0.4 – 40 mAs
60 – 68 kVp	0.4 – 40 mAs
70 – 78 kVp	0.4 – 40 mAs
80 – 88 kVp	0.4 – 40 mAs
90 – 100 kVp	0.4 – 40 mAs

**Every system includes:** AE-1020HF High Frequency Portable X-Ray, power cord, exposure switch, carrying case, operators manual.

---

## Contact Information:



**apexx Veterinary Equipment**  
**14 Inverness Drive. East, Suite C-124**  
**Englewood, CO 80112**

**Ph: 800.211.0036**  
**Email: [help@apexx360.com](mailto:help@apexx360.com)**

**[www.apexx360.com](http://www.apexx360.com)**

\* Information and specifications subject to change without notice.